


# Coronis 3MP (MFGD-3621)

3 MegaPixel premium diagnostic display system



Coronis 3MP features a 20.8" grayscale flat panel with a resolution of 1536 x 2048 for perfect image representation. With a contrast ratio of 600:1 and high luminance, Coronis 3MP offers medical professionals a highly accurate solution for a host of medical imaging applications including PACS, computed tomography, MRI, angiography and chest radiology. The 3MP display system can be used both in dark reading rooms as well as in ambient light environments. The integrated Coronis Portrait Accelerator assures full performance in portrait mode.

**BARCO**

Visibly yours

# Technische Daten

<b>Zulassungen</b>	<ul style="list-style-type: none"><li>• CE, UL60950, CSA C22.2 No 60950 (c-UL), CB , IEC 60950, DEMKO EN 60950</li></ul>
<b>Unterstützte Auflösungen</b>	<ul style="list-style-type: none"><li>• 2048 x 1536 landscape</li><li>• 1536 x 2048 portrait (features Portrait Accelerator)</li><li>• VGA boot (1 to 1 pixel mapping)</li></ul>
<b>Display-Technologie</b>	<ul style="list-style-type: none"><li>• TFT AMLCD dual domain IPS</li><li>• Grayscale, 3 subpixels per pixel</li></ul>
<b>Ergonomie</b>	<ul style="list-style-type: none"><li>• Tilt &amp; swivel base with P/L rotation</li><li>• 100 mm VESA mounting standard</li><li>• AR coated front cover</li></ul>
<b>Auflösung</b>	3MP (2048 x 1536 / 1536 x 2048)
<b>Leuchtdichte</b>	500 cd/m <sup>2</sup> (145.9 fL) calibrated, 700 cd/m <sup>2</sup> (204 fL) max.
<b>Dunkelraumkontrast</b>	600:1 typical
<b>DVI</b>	Dual DVI (complying to VESA DVI)
<b>Stabilisierung, Kalibrierung und automatisierte QS</b>	Embedded I-Guard
<b>Analog</b>	No
<b>USB</b>	USB-hub with 1 up- and 2 downstream ports